

CellCover

Basic Protocol No. 3

WWW.ANACYTE.COM

ANACYTE^{••}
LABORATORIES

CellCover

Basic protocol No. 3

CellCover exerts its stabilizing effect very fast. DNA, RNA, protein stay in place in a close to native condition, without crosslinking activity! Harsh treatments with alcohol, acetone or even formaldehyde can be avoided.

Basic protocol for tissue specimens:

- ① **Cut tissue into approximately 5 mm pieces**
- ② **CellCover does not penetrate tight junctions, thus encapsulated organs must be open across the largest diameter**
- ③ **Place tissue in CellCover (at least 10x volumes) and store at 4°C until use**
- ④ **Change CellCover at least once after 4-24 hours**
- ⑤ **Proceed according to experimental design e.g. RNA isolation**

Possible downstream applications:

- **Laser capture microscopy**
- **ISH: in situ hybridization (RNA as well as DNA FISH and CISH!)**
- **Batch/ single cell transcriptome analysis**

For questions concerning experimental strategies and special applications of Anacyte's products, please contact our support:

TEL +49 (0)40 60941702-0

E-MAIL support@anacyte.com

WEB www.anacyte.com

Anacyte Laboratories GmbH

Saseler Bogen 3

22393 Hamburg, Germany

ANACYTE
LABORATORIES